01 | 2019 3. JAHRGANG

THERAPY

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LITERATUR

- [1] **Aman J. E.** (2015) The effectiveness of proprioceptive training for improving motor function: a systematic review, *Front Hum Neurosci.* 8:1075.
- [2] **An M.** (2011) The effects of exercise-based rehabilitation on balance and gait for stroke patients: a systematic review. *J Neurosci Nurs.* 43(6):298-307.
- [3] **Batchelor F. et al.** (2012) Falls after stroke. *International Journal of Stroke.*
- [4] **Bonan I. V.** (2004) Reliance on visual information after stroke. Part II: Effectiveness of a balance rehabilitation program with visual cue deprivation after stroke: a randomized controlled trial. *Arch Phys Med Rehabil*;85:274-8.
- [5] **Cheng F.-Y.,** (2014) Factors Influencing Turning and Its Relationship with Falls in Individuals with Parkinson's Disease. *PLoS ONE* 9(4): e93572.
- [6] **Davies B. L.** (2015) Neurorehabilitation Strategies Focusing on Ankle Control Improve Mobility and Posture in Persons With Multiple Sclerosis. *J Neurol Phys Ther.* 39(4):225-32.
- [7] **de Haart M.** (2004) Recovery of standing balance in postacute stroke patients: a rehabilitation cohort study. *Arch Phys Med Rehabil* 85:886-95.
- [8] **Dohle Ch.** (2015) Rehabilitation der Mobilität nach Schlaganfall (ReMoS). S2e-Leitlinie. *Neurologie & Rehabilitation* 7.
- [9] **EBRSR** (2016) The Stroke Rehabilitation Clinician's Handbook. Download am 07.05.2018 unter <http://www.ebrsr.com/clinician-handbook>.
- [10] **Gandolfi M.** (2015) Sensory integration balance training in patients with multiple sclerosis: A randomized, controlled trial. *Mult Scler.* 21(11):1453-62.
- [11] **Geurts A. C.** (2005) A review of standing balance recovery from stroke. *Gait Posture* 22(3):267-81.
- [12] **Horak F. B.** (2006) Postural orientation and equilibrium: what do we need to know about neural control of balance to prevent falls. *Age Ageing.* 35 Suppl 2:ii7-ii11.
- [13] **Huber M.** (2016) „Posturale Kontrolle – Grundlagen.“ *Neuroreha* (Thieme-Verlag). 8: 158-162.
- [14] **Huber M.** (2018) Balancepad – wissen wir wie's wirkt? *physio-praxis* 2018; 16(5):30-31.
- [15] **Jöbges M.** (2004) Repetitive training of compensatory steps: a therapeutic approach for postural instability in Parkinson's disease *J Neurol Neurosurg Psychiatry.*75:1682-1687.
- [16] **Kitatani R. et al.** (2016) Ankle muscle coactivation during gait is decreased immediately after anterior weight shift practice in adults after stroke. *Gait & Posture* 45:35-40.
- [17] **KNGF** (2014) Clinical Practice Guideline for Physical Therapy in patients with stroke. <https://www.fysionet-evidencebased.nl/index.php/kngf-guidelines-in-english>.
- [18] **Mansfield A.** (2011) Training rapid stepping responses in an individual with stroke. *Phys Ther.* 91:958-969.
- [19] **Marigold D. C.** (2005) Exercise Leads to Faster Postural Reflexes, Improved Balance and Mobility, and Fewer Falls in Older Persons with Chronic Stroke. *JAGS* 53:416-423.
- [20] **Matsuda P.** (2011) Falls in Multiple Sclerosis, *PM&R* 3;(7) : 624-632.
- [21] **Robinovitch S. N.** (2013) Video capture of the circumstances of falls in elderly people residing in long term care: an observational study. *Lancet* 381(9860): 47-54.
- [22] **Shumway-Cook A., Woollacott M. J.** (2016) *Motor Control* 5th Edition Lippincott Raven.
- [23] **Sitsupadol P. et al.** (2006) Training of balance under single and dual-task conditions in older adults with balance impairment. *PhysTher.* 86:269-281.
- [24] **Smania N.** (2008) Rehabilitation of sensorimotor integration deficits in balance impairment of patients with stroke hemiparesis: a before/after pilot study. *Neurol Sci.* 29:313-319.
- [25] **Sparrow D.** (2016) Highly Challenging Balance Program Reduces Fall Rate in Parkinson Disease. *J Neurol Phys Ther.*40(1):24-30.
- [26] **Stroke Foundation** (2017) *Clinical Guidelines for Stroke Management.* Melbourne Australia. <https://informme.org.au/Guidelines/Clinical-Guidelines-for-Stroke-Management-2017>.
- [27] **Tasseel-Ponche S.** (2015) Motor strategies of postural control after hemispheric stroke. *Neurophysiol Clin.* 45(4-5):327-33.
- [28] **van Duijnhoven H. J.** (2016) Effects of Exercise Therapy on Balance Capacity in Chronic Stroke: Systematic Review and Meta-Analysis. *Stroke.* 47(10):2603-10.
- [29] **Veerbeek J. M.** (2014) What Is the Evidence for Physical Therapy Poststroke? A Systematic Review and Meta-Analysis. *PLoS ONE* 9(2): e87987.

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LITERATUR

[1] Ada, L./Canning, C./Low, S.-L.: Stroke patients have selective muscle weakness in shortened range, Brain 2003, 724-731.

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Effects of a device-based over-ground gait training

Regaining ability to walk again is one of the most important goals in the field of neurological rehabilitation. In this paper, we describe the use of the motorized assistive gait trainer device THERA-Trainer e-go and the corresponding effects of an over-ground gait training in neurological patients. Results suggest that the proposed training paradigm may successfully ascertain training conditions that correspond well to patient capabilities and therefore may potentially be useful in improving over-ground walking during gait rehabilitation of neurological individuals.

Jakob Tiebel

Background

Regaining ability to walk again is one of the most important goals in the field of neurological rehabilitation.

In ambulatory patients with neurological disease often an impaired coordination of walking and a significant reduced walking speed and endurance can be observed [2]. Only about 7% of individuals discharged from inpatient rehabilitation are able to achieve average walking speed and distance, required for unrestricted mobility in the community [7].

In clinical practice training of walking on a task specific and repetitive basis, seems to be most effective to improve patient's gait performance [2, 3].

However, next to specific interventions, like speed-dependent treadmill training, walking should be trained under conditions of everyday life [4, 5, 7].

In order to practice an intensive training of over-ground walking under fall-safe conditions, a mobile over-ground gait trainer system that supports a patient much as needed can be a helpful and decisive advantage [1, 8].

Objective

Main objective of this case series study was to report feasibility and observed effects of a device-based over-ground walking training with the gait trainer device THERA-Trainer e-go.

Patients

The study was conducted at the Stroke Rehabilitation Hospital RehaNova in Cologne, Germany. Three patients with neurological disease were recruited by physiotherapists according to inclusion criteria between March 2014 and April 2014. Generally, we included inpatient adults with an upper motor-neuron syndrome (UMNS), maximal 12 weeks post-injury and a mobility impairment determined by using the Functional Ambulation Categories (FAC). At study onset patients at least had to be able to walk minimum 10 meters with assistance of one helping person, which corresponds to an FAC score ≥ 2 .

In the following we report on the patients with average post-injury time of 8 weeks. All three

patients met inclusion criteria and agreed to participate before intervention started.

Subject 1: a 66-year-old man with post-stroke right hemiparesis, 5 weeks post-injury and FAC 2. Subject 2: a 70-year-old man with Guillain-Barre syndrome, 12 weeks after hospital admission and FAC 3.

Subject 3: a 57-year-old woman with post-stroke left hemiparesis, 8 weeks post-injury and FAC 2.

Intervention

Patients trained walking for three weeks with the THERA-Trainer e-go device, 3 - 5 times per week in addition to conventional physiotherapy.

Standard physiotherapy training was 45 min/day, followed by the training with the gait trainer device for a duration according to each patient's endurance. The routine motor rehabilitation used a task-oriented and goal-oriented approach that included gait-preparatory interventions and balance exercises while sitting and standing. The treatment was modified according to the patient's needs and the rehabilitation teams' goals. Walking aids and orthoses have been prescribed according to individual needs of the patients.

Description of the device

The THERA-Trainer e-go is a motorized assistive gait trainer device developed by medica Medizintechnik GmbH (www.thera-trainer.com).

Two independent motorized wheels, one on each side of the base platform of the device, provide necessary actuations to exercise the essential movement manoeuvres (moving forward, accelerating/decelerating and turning left/right) with adjustable linear and angular velocity. On the same base platform helical springs, one on each side, are housed in steel cylinders. Between the helical spring and the walls of the steel cylinder, there is a resistance adjustment ring with a handle. By displacing the adjustment ring upwards, the bending length of the spring becomes shorter, thus making the spring stiffer. Conversely, by lowering the adjustment ring the bending length of the spring is longer, thus making the spring more compliant. Each spring connects to a vertical rod that is coupled with supporting horizontal bars and the table


Fig 2. Schematic representation of Gait Balance Trainer. The helical springs allow the standing frame to comply in sagittal and frontal plane while interacting with the subject during training which assures natural pelvis movement [6].

top to constitute a supporting frame. During training the individual is embraced around pelvis with harness and attached to the supporting frame on each side in a way that minimizes relative movement between the supporting frame and the individual, but at the same time to allow natural pelvis movement.

Therefore, if the individual leans forward or backward and/or sideways, the supporting frame follows the movement, which provokes the helical spring to bend, thus providing horizontal support. The level of supporting force is selected such that it provides minimal needed assistance.


The speed and direction of the e-go is therapist controlled with a manual control module. The device measures time and walking distance of each training session. Maximal walking speed is 3 km/h. The device does not offer body weight support; however, leather harness suspended from horizontal rods can easily accept whole weight of an individual in case of sudden vertical drop due to insufficient strength of leg muscles, to ensure fall-safe training conditions.


Measurements

The following tests were performed: walking capacity (measured with 2-Minute-Walk-Test, 2-MWT),

walking speed (measured with 10-Meter-Walk-Test, 10-MWT) and Functional Ambulation Category (FAC), to evaluate ambulation ability. The severity of activity limitation in ADL was evaluated using the Functional Independence Measure.

Measurements were taken through within-session and across-session changes at baseline and after the intervention period of 3 weeks with a pre-test (before training session), post-test (after training session) and retention-test (24h after training session). In addition, we also measured


the walking distance (average walking distance of the training sessions) and the walking speed (average walking speed of the training sessions).

Results

No adverse events occurred during training. Physiotherapists did not report any significant problems with the training device handling, and there was no additional staff required to set up the participant in the device.

As expected, all patients improved from baseline to end of study measurement in walking speed, walking distance, FAC and FIM score.



Subject 1: After three weeks of training, the patient's walking ability changed from FAC 2 to 4. While at study onset the patient was able to walk


with help of one therapist, after three weeks of training the patient was able to walk independently with a crane inside the hospital. At study end, there was a significant improvement of 31.4% (FIM: 45 to 79) in ADL independency.

Walking speed, increased about 0.46 m/s (wk3: 0.82 m/s - wk1: 0.36 m/s) in 10-MWT and walking capacity increased about 50.7 m (wk3: 89.1 m - wk1: 38.4 m) measured with 2-MWT. Comparison based on the first with last week's retention test results.


The evaluation of the averaged training results shows a continuous extension of walking distance and walking speed over the course of three weeks training. From first to second week walking speed increased about 0.05 m/s (wk2: 0.38 m/s - wk1: 0.33 m/s) and from second to third week about 0.2 m/s (wk3: 0.58 m/s - wk2: 0.38 m/s).

Walking distance increased from first to second week about 142 m (wk2: 278 m - wk1: 136 m) and from second to third week about 226 m (wk3: 504 m - wk2: 278 m). The overall changes comparing the results from first and third week, were 0.25 m/s (wk3: 0.58 m/s - wk1: 0.33 m/s) walking speed and 368 m (wk3: 504 m - wk1: 136 m) walking distance.

Subject 2: After three weeks of training, the patient's walking ability changed from FAC 3 to 4. While at study onset the patient was able to walk supervised by a therapist, after three weeks of training the patient was able to walk independently with a cane inside the hospital. At study end, there was a significant improvement of 22.2% (FIM: 71 to 95) in ADL independency.

Walking speed, increased about 0.23 m/s (wk3: 0.81 m/s - wk1: 0.58 m/s) in 10-MWT and walking capacity increased about 31.4 m (wk3: 89.1 m - wk1: 57.7 m), measured with 2-MWT. Results based on comparison of the first with last week's retention test results.

The evaluation of the averaged training results shows a continuous extension of walking distance and walking speed over the course of three weeks training. From first to second week walking speed increased about 0.15 m/s (wk2: 0.71 m/s - wk1: 0.56 m/s) and from second to third week about 0.04 m/s (wk3: 0.75 m/s - wk2: 0.71 m/s).

Walking distance increased from first to second week about 412 m (wk2: 767 m - wk1: 355 m) and from second to third week about 276 m (wk3: 1043 m - wk2: 767 m).

The overall changes comparing the results from first and third week, were 0.19 m/s (wk3: 0.75 m/s - wk1: 0.56 m/s) walking speed and 688 m (wk3: 1043 m - wk1: 355 m) walking distance.

Subject 3: After 3 weeks of training, the patient's walking ability changed from FAC 2 to 4. While at study onset the patient was able to walk with help of one therapist, after three weeks of training the patient was able to walk independently on even surfaces with a walker inside the hospital.

The evaluation shows a continuous extension of walking distance and walking speed.

Walking speed, increased about 0.30 m/s (wk3: 0.76 m/s - wk1: 0.46 m/s) in 10-MWT and walking capacity increased about 35.7 m (wk3: 85.9 m - wk1: 50.2 m), measured with 2-MWT. Results based on comparison of the first with last week's retention test results.

The evaluation of the averaged training results shows a continuous extension of walking distance and walking speed over the course of three weeks training. From first to second week walking speed increased about 0.18 m/s (wk2: 0.56 m/s - wk1: 0.38 m/s) and from second to third week about 0.13 m/s (wk3: 0.69 m/s - wk2: 0.56 m/s).

Walking distance increased from first to second week about 188 m (wk2: 383 m - wk1: 195 m) and from second to third week about 436 m (wk3: 819 m - wk2: 383 m).

The overall changes comparing the results from first and third week, were 0.31 m/s (wk3: 0.69 m/s - wk1: 0.38 m/s) walking speed and 624 m (wk3: 819 m - wk1: 195 m) walking distance.

Discussion

In this series case study, we described the use of the motorized assistive gait trainer device THERA-Trainer e-go and the corresponding effects of an over-ground walking training in neurological in-patients.

Results suggest that the proposed training paradigm may successfully ascertain training conditions that correspond well to patient capabilities and therefore may potentially be useful in improving over-ground walking during gait rehabilitation of neurological individuals.

We demonstrate the feasibility and advantage of using the within- and across-session changes for evaluating the improvements during clinical gait rehabilitation with the gait trainer device.

In conclusion, the THERA-Trainer e-go appears to be a comfortable and easy to use gait trainer device, that confers a sense of freedom and safety to subjects during over-ground gait training. That helps to increase the intensity of training. During over ground gait training the device can safely deliver assistance to patients with risk of falling and helps to improve ability to walk independently. That also seems to have an impact on patients' independence in activities of daily living.

Acknowledgements

We thank all patients and their caregivers who participated in this study. We gratefully acknowledge the support of the clinical staff at the RehaNova and all those involved in this project.

Consent for publication

Subjects of the study gave consent to use and publish data in such way that anonymity will be assured.

LITERATURE

[1] **Bizovičar N, Matjačić Z, Stanonik I, Goljar N** (2017). Over-ground gait training using a motorized assistive device in patients with severe disabilities after stroke. *International Journal of Rehabilitation Research*, 40:46-52.

[2] **Hollands KL, Pelton TA, Tyson SF, Hollands MA, van Vliet PM** (2012). Interventions for coordination of walking following stroke: systematic review. *Gait Posture* 35:349-359.

[3] **Lohse KR, Lang CE, Boyd LA** (2015). Is more better? Using meta-data to explore dose-response relationships in stroke rehabilitation. *Stroke* 45:2053-2058.

[4] **Mehrholz J** (2016) *Neurorehabilitation von Stand und Gang*. [Neurorehabilitation of standing and walking.] In Platz, Th. Update Neurorehabilitation. Bad Honnef: Hippocampus Verlag.

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[6] **Olensek A, Oblak J, Cikajlo I, Novak P, Jere K, Matjacic Z** (2012). Adaptive dynamic balance training during overground walking with assistive device. *Proceedings of the IEEE RAS and EMBS International Conference on Biomedical Robotics and Biomechanics*, (July 2015), 1066-1070.

[7] **Phan PL, Blennerhassett JM, Lythgo N, Dite W, Morris ME** (2013). Over-ground walking on level and sloped surfaces in people with stroke compared to healthy matched adults. *Disabil Rehabil* 35:1302-1307.

[8] **Tiebel J** (2018). Modern gait rehabilitation. Putting hybrid gait trainers to the test. *THERAPY Magazine* 1:30-37.

How to cite this article:

Tiebel, J (2018). Effects of a device-based over-ground walking training with a novel motorized assistive gait trainer device THERA-Trainer e-go. A case series study.

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LITERATUR

- [1] Duncan, P. W./Studenski, S./Chancler, J./Prescott B.: Functional reach: Predictive validity in a sample of elderly male veterans, in: Journal of Gerontology 47, 1992, 3, 93-98.
- [2] Duncan, P. W./Weiner, D. K./Chancler, J./Studenski, S.: Functional reach: A new clinical measure of balance, in: Journal of Gerontology 45, 1990, 6, 192-197.
- [3] Tielbel J. (2018): Intensive continuous therapy for paraplegia, THERAPY Magazine, 2(1), 38-41.
- [4] Tielbel, J. (2018): Back to everyday life, step by step, THERAPY Magazine, 2(1), 42-43.
- [5] Tielbel, J. (2018): Practice walking by walking, THERAPY Magazine, 2(1), 20-24.
- [6] Tielbel, J. (2018): Putting hybrid gait trainers to the test, THERAPY Magazine, 2(1), 30-37.
- [7] Yardley, L./Beyer, N./Hauer, K./Kempen, G./Piot-Ziegler, C./Todd, C. (2005): Development and initial validation of the Falls Efficacy Scale International (FES-I). Age and Ageing, 34(6), 614-619.

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LITERATUR

Santos L, et al. Effects of progressive resistance exercise in akinetic-rigid Parkinson's disease patients: a randomized controlled trial. Eur J Phys Rehabil Med 2017 Mar 13. DOI: 10.23736/S1973-9087.17.04572-5

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
The objective of the ReMoS guideline was the investigation of effects of rehabilitation with the goal of improving the function of the lower extremity. Clinical experience has shown that therapeutic interventions are not equally effective at all times post stroke.

Recommendations from the German Guideline for „Rehabilitation of Mobility after Stroke”

Guidelines reflect the state of the art and generate a high degree of confidence that their recommendations support optimal treatment. Since 2015 the German Neurorehabilitation Society (DGNR) provides an evidence-based guideline for “Rehabilitation of Mobility after Stroke” (ReMoS). Unfortunately, until now it is available exclusively in German. To overcome the language barrier this article covers main recommendations of the guideline in English.

Jakob Tiebel

Background and key issues

The ReMoS guideline was jointly developed by neurologists and physiotherapists working in the field of neurorehabilitation. It is the first stroke rehabilitation guideline that is strictly oriented on target criteria and in which authors make differentiated recommendations for patients in the sub-acute or chronic stage after stroke and for initially non-ambulatory and ambulatory patients (ReMoS Working Group, 2015 at 355).

The guideline focuses on the following key issues:

1. What rehabilitation measures for patients with a stroke or hemiparesis after stroke have been proven to lead to an improvement in
 - a) walking ability
 - b) risk of falling and balance
 - c) walking speed and/or walking distance?

2. What rehabilitation methods can be recommended for improving the target criteria mentioned above?


In the developing process also other national and international guidelines in the field of stroke rehabilitation have been taken into account, provided that these guidelines were devoted to the same interventions. But a direct adoption of recommendations was not possible, since ReMoS gives recommendations for interventions with regard to specific target criteria, whereas in the other guidelines reviewed recommendations are not differentiated according to that. Thus, the focus on target criteria is an important differentiator compared to other stroke guidelines in the world and shall make it much easier to use it in clinical practice (ReMoS Working Group, 2015 at 355-57).

Classification of different stages post stroke

Clinical experience has shown that therapeutic interventions are not equally effective at all times post stroke. In order to better assess the differential efficacy of interventions, the authors investigated their efficacy in the acute, subacute or chronic stages after stroke.

However, the acute, subacute and chronic phases are defined differently in the field of neurological rehabilitation. Time windows normally range between 3 weeks and 12 months post stroke.

For the purposes of the ReMoS guideline the researchers defined the time-windows of 3 weeks and 6 months as boundaries between the acute, sub-

The focus on target criteria is an important differentiator compared to other stroke guidelines and shall make it easier to use it in clinical practice.

acute and chronic phase of rehabilitation (ReMoS Working Group, 2015 at 356-57).

General approach

The following approach was adopted for the development of the ReMoS guideline (ReMoS Working Group, 2015 at 358):

1. Systematic literature search using defined criteria
2. Evaluation of individual papers
3. Summary of individual papers relating to a question and evaluation of the baseline quality of the evidence
4. Critical summary while considering the appropriateness and clarity of the evidence, resulting in the final quality of the evidence (GRADE scheme)
5. Practical consideration (e.g. risk of the therapy procedure) and development of a recommendation from this.

The systematic search was based on the following PICO key question (ReMoS Working Group, 2015 at 358):

Patients:

In-patients with a stroke or a hemiparesis after stroke...

Intervention:

... specific rehabilitation therapy (physiotherapy, occupational therapy, acupuncture, electrical stimulation, robot-assisted gait therapy, treadmill training, biofeedback therapy, implantation of medical devices, botulinum toxin injection, orthotic supplies and other rehabilitation therapy)...

Comparison:

... with varying doses or different contents...

Outcomes:

... lead to improvement in strength or mobility in the lower extremities or improvement in balance, walking, mobility or reduction in number of falls?


Literature search

The literature search was done with the program Reference Manager in the Medline (Pubmed) and PEDro databases. The last systematic search took place in June 2012. The PICO search algorithm included the same terms for all specified databases, but a syntactic change of the search terms was necessary depending on the version of the internet database. Particularly in the case of Medline (Pubmed) the authors stated, that during the process of literature search the internal search algorithm changed considerably a number of times. Therefore, before each literature search they ran several test searches with an orientation triage of the results and an adjustment of the search algorithm (ReMoS Working Group, 2015 at 358).

Interventions included

The researchers carried out an open literature search with regard to interventions that selected references according to target criteria and patient cohort only. As a result, also studies were found in which one of the target criterion was not primarily intended and a plausible, specific effect mechanism was not directly established (ReMoS Working Group, 2015 at 359).

However, the objective of the ReMoS guideline was the investigation of effects of rehabilitation (such as physiotherapy, occupational therapy, acupuncture, electrical stimulation, assisted gait therapy, biofeedback therapy and medication) with the goal of improving the function of the lower extremity. Tests were being carried out to determine whether there was at least one study in which this intervention was investigated as a primary intervention. In case of this, the effect of the intervention was also considered in studies where this intervention was treated as a control intervention. Studies that included purely acute medical aspects therapy were not taken into account (ReMoS Working Group, 2015 at 359).

Studies included

The researchers included only randomized controlled trials, randomized crossover studies, systematic reviews and meta-analyses. In case of inclusion of cross-over studies always only the first phase of the intervention was considered. And exclusively articles published in English or German on adult stroke patients were included (ReMoS Working Group, 2015 at 359).

Additionally, sample sizes were decisive for the inclusion of a study in the guideline, as small case numbers could easily lead to an inappropriate distortion of the results. Therefore, randomized trials with less than 10 included patients per group were not considered. Studies with fewer than 20 patients included per group were generally downgraded in the quality ranking even if the quality was otherwise good (ReMoS Working Group, 2015 at 359). Nevertheless, as not enough studies with the desired quality were available for all interventions, the researchers made exceptions to this rule. If there were several "small" studies with less than 10 patients per group, but a meta-analysis with more than 10 patients per group, the meta-analysis was evaluated. If there were no "larger" studies with more than 10 patients per group for an interven-

tion available, but one or more smaller studies with less than 10 patients per group, the "smaller" studies were considered in the work of the guideline (ReMoS Working Group, 2015 at 359).

Target criteria


The classification of the target criteria was limited to the following clinically most relevant points by Consensus of the guideline working group (ReMoS Working Group, 2015 at 360):

- Walking ability
- Walking speed
- Walking endurance (in relation to distance)
- Risk of falling (Balance)

This classification forms the structure of the guideline according to these target criteria. General motor scales such as the Fugl-Meyer Assessment, Rivermead Motor Assessment or total ADL scales such as Barthel Index and FIM are disregarded explicitly. Papers whose endpoints were exclusively physiological parameters, e.g. muscle strength or kinematic parameters, were also not included in the guideline (ReMoS Working Group, 2015 at 360-62).

Number of studies included

The Researchers included a total of 199 papers, of which 188 were clinical trial papers and 11 were systematic reviews (ReMoS Working Group, 2015 at 362).

Evaluation of individual papers

All studies were initially evaluated irrespective of the guideline question.

By rating the methodology the authors determined the validity of studies and classified them according to the Oxford classification, (Version 1 of the CEBM classification of March 2009) (ReMoS Working Group, 2015 at 363-64).

In order to produce a practical recommendation, studies were allocated as far as possible to overarching therapeutic approaches, for example (ReMoS Working Group, 2015 at 364):

- Conventional gait training
- Treadmill training with patients in the subacute stage
- Treadmill training with patients in the chronic stage
- End-effector-based devices (e.g. Gait Trainer)
- Exoskeleton-assisted devices (e.g. Lokomat)
- Training for improved strength and endurance ergometer training
- Specific balance training
- Sensory feedforward stimulation without feedback (except electrical stimulation)
- Cognitive therapy strategies
- Comprehensive physiotherapy treatment concepts
- Electrical stimulation
- Aids
- Central stimulation
- Drug therapy
- Positioning

This type of structuring directly affected evidence assessment and recommendations. In the case of different device-based procedures, for example, the researchers did not look at every industry model in isolation. They examined the overall principle of treatment (e.g. end-effector versus exoskeleton). It should be remembered that this type of classification is subject to a certain degree of subjectivity. But the structure makes much sense and helps therapists to better find recommendations for different intervention types. The selected classifications were agreed in the whole guideline working group (ReMoS Working Group, 2015 at 364).

In addition, therapy organization concepts were evaluated as follows (ReMoS Working Group, 2015 at 364):

- Organization of the therapy sessions
- Training intensity, additional physiotherapy
- Training organization
- Circuit class training

- Cross-organizational concepts
- Early supported discharge with multidisciplinary rehabilitation in the domestic environment
- Organisational concepts in the domestic environment based devices (e.g. gait trainers)

This structure also enabled a comprehensive discussion of possible principles of treatment.

In the summary evaluation, the formal baseline quality of the evidence for an intervention has been determined in the four-stage GRADE classification, starting with the methodological quality of the individual papers (ReMoS Working Group, 2015 at 365).

Table 1: GRADE System – Strength of Recommendation (Guyatt et al., 2008, and GRADE Working Group)

Quality of Evidence	Definitions
High	Further research is very unlikely to change our confidence in the estimate of effect.
Moderate	Further research is likely to have an important impact on our confidence in the estimate of effect.
Low	Further research is very likely to have an important impact on our confidence in the estimate of effect.
Very low	Any estimate of effect is uncertain.


The recommendations are formulated on the basis of the evaluation of the quality of the evidence in the consensus procedure while taking into account potential risks and clinical benefits. Similar to other DGNR guidelines (for example the guideline for the rehabilitation of the upper extremities post stroke from Platz et al.) the following recommendation levels are used, which can be formulated both positively and negatively (ReMoS Working Group, 2015 at 365):

Table 2: Recommendations

Recommendation level	Positive formulation	Negative formulation
A	ought	ought not
B	should	should not
C	can be carried out	"cannot be recommended"
No recommendation		

In the formulation of the open recommendation, a distinction was made as to whether the recommendation takes effect on the basis of insufficient strength of evidence, for example a very small study, or due to lack of effects, including larger studies. Only in the first case is a "can" recommendation (0) expressed. In the second case, always a "no recommendation" was expressed by the authors

All evaluations of evidence and the derivation of recommendations were done separately for the interventions. That made possible to summarize recommendations for the individual target criteria in different chapters.

without specifying a recommendation grade (ReMoS Working Group, 2015 at 365).

Summary evaluation of interventions

As shown above, all evaluations of evidence and the derivation of recommendations were done separately for the interventions. That made possible to summarize recommendations for the individual target criteria in different chapters. The following overview should enable the selection of suitable interventions in a concrete case under a specific question (ReMoS Working Group, 2015 at 366).

Walking ability

A high number of steps carried out (number of repetitions) is a significant principle of treatment to improve walking ability. When selecting the therapies to use, a decision must be made when therapy is initiated as to whether the patients are able to walk or not. In the case of non-ambulatory patients, the following recommendations can be expressed (ReMoS Working Group, 2015 at 481):

Table 3: Attainment of walking ability in non-ambulatory patients (ReMoS Working Group, 2015 at 482)

	Subacute	Chronic
B (should):	– Intensive gait training, if available and feasible, including end-effector-based devices (e.g. Gait Trainer)	
0 (can):	– Intensive gait training, if available and feasible, including treadmill or exoskeleton-assisted devices (e.g. Lokomat) – Cyclical multichannel stimulation to generate gait-similar leg movements in the paretic leg when lying – Additional electroacupuncture – For patients with neglect : Specific neglect training	

Especially in the subacute stage, the high number of repetitions achieved, if available, including a training on an end-effector device (e.g. gait trainer). Alternatively, gait training can also be performed with an exoskeleton (e.g. Lokomat) or treadmill (ReMoS Working Group, 2015 at 481). Compared to the treatment of non-ambulatory patients, the

device-support can no longer be seen as as preferential for at least assisted ambulatory patients (ReMoS Working Group, 2015 at 481):

Table 4: Improvement in walking ability in (restricted) ambulatory patients (ReMoS Working Group, 2015 at 482).

	Subacute	Chronic
B (should):	– Intensive gait training: conventional or including treadmill (possibly progressive)	– For patients with spastic equinovarus deformity: injection of botulinum toxin to reduce the use of aids
0 (can):	– Task-related training with movement imagery – Use of walking aids – Combination therapy of end-effector-based gait training (e.g. Gait Trainer) with functional electrical stimulation – Needle acupuncture including electro-acupuncture during intensive rehabilitation	– Support of treadmill training with VR
– B (should not):		– For patients with spastic equinovarus deformity: Thermo-coagulation of the tibial nerve

In these patients, conventional or device supported training (in this case usually with the treadmill) can be made relatively freely dependent on the possibilities of the facility, without significant differences in the result. If the treadmill is used, it should be used as progressively as possible (increasing the speed) (ReMoS Working Group, 2015 at 481).

Walking speed

If the aim is to improve the walking speed in ambulatory patients, progressive aerobic treadmill training can be considered as a gold standard. Also intensive walking practice with and without treadmill or an intensive supervised home exercise program can be strongly recommended. The same also applies to gait training with acoustic feedback (ReMoS Working Group, 2015 at 482).

For patients with initial (limited) walking ability, a step-synchronous stimulation of flexoreflex- afferent nerve pathways during walking training can be recommended in order to achieve a higher walking speed (ReMoS Working Group, 2015 at 482).

However, the use of a gait trainer or acoustic training in severely affected patients (ReMoS Working Group, 2015 at 482).
stimulation can also improve the effects of gait Working Group, 2015 at 482).

Table 5: Improvement in speed in (restricted) ambulatory patients (ReMoS Working Group, 2015 at 482)

	Subacute	Chronic
A (ought):	<ul style="list-style-type: none"> - Task-related progressive endurance training (in implementation treadmill or progressive circuit class training) 	
B (should):	<ul style="list-style-type: none"> - Intensive gait training without treadmill or - Intensive gait training including treadmill or - Intensive supervised home exercise programme (strengthening, endurance, balance) with progression - Gait training with stimulation of flexor reflex afferents 	<ul style="list-style-type: none"> - Orthosis with electrical stimulation of the peroneal nerve (indirect effect)
0 (can):	<ul style="list-style-type: none"> - Intensive progressive task-related training - Task-related training with movement imagery - Walking training with end-effector or exoskeleton gait trainer if device is available - Isokinetic strength training - Strength/endurance training - Gait training with acoustic stimulation - Acoustic feedback while walking - Feedback/reinforcement (daily time measurement while walking with reinforcing feedback) - Combination therapy of end-effector-based gait training (e.g. Gait Trainer) with functional electrical stimulation - Electroacupuncture - Ankle brace - Early orthopaedic shoe - In case of severe arm paresis: arm sling 	<ul style="list-style-type: none"> - Intensive progressive task-related training - Task-related endurance training, e.g. progressive aerobic treadmill training - Combination of treadmill training and variable gait training on the ground - Gait training with end-effector-based gait trainer (e.g. Gait Trainer) if device is available - Additional backward walking training - Task-related strength training - Additional balance training - Strength feedback training - Additional VR-based training - Additional dual-task training - Movement observation - TENS on tendon transition of the spastic medial gastrocnemius - TENS on acupuncture points before task-oriented training - Ankle brace - Orthosis with electrical stimulation of the peroneal nerve (direct effect) - Toe separator with/without shoe - Repetitive magnetic stimulation in combination with task-oriented training
- B (should not):		<ul style="list-style-type: none"> - For patients with spastic equinovarus deformity: Thermocoagulation of the tibial nerve

Walking Distance

The following recommendations can be expressed to improve distance:

Table 6: Improvement in distance in (restricted) ambulatory patients (ReMoS Working Group, 2015 at 482).

	Subacute	Chronic
A (ought):	<ul style="list-style-type: none"> - Task-related progressive endurance training 	
B (should):	<ul style="list-style-type: none"> - Intensive supervised home exercise programme (strengthening, endurance, balance) with progression - Task-related aerobic treadmill training - Intensive gait training including treadmill 	<ul style="list-style-type: none"> - Task-related endurance training, e.g. progressive aerobic treadmill training - Orthosis with electrical stimulation (indirect effect)
0 (can):	<ul style="list-style-type: none"> - Task-related training with movement imagery - Gait training with end-effector or exoskeleton gait trainer if device is available - Strength/endurance training - Additional strength/endurance training - Multichannel functional electrical stimulation during gait training 	<ul style="list-style-type: none"> - Intensive progressive task-related training - Combination of treadmill training and variable gait training on the ground - Gait training with end-effector gait trainer if device is available - Ergometer training - Task-related strength training - Strength feedback training - Movement observation - Orthosis with electrical stimulation (direct effect)
- B (should not):		<ul style="list-style-type: none"> - Functional electrical stimulation with percutaneous wire electrodes

Improving the walking distance in walking patients requires in particular an increase in cardiopulmonary performance. The gold standard for this is aerobic endurance training, which should be carried out in relation to the task at hand, e.g. as progressive aerobic treadmill training (ReMoS Working Group, 2015 at 483).

Here, as well, an intensive supervised home exercise program can be highly recommended, as well as an intensive progressive task-oriented training or an isolated strength endurance training. Also in this case additional effects can be achieved through electrostimulation procedures, especially for severely affected patients (ReMoS Working Group, 2015 at 483).

Risk of falling, balance

The following recommendations can be expressed to improve balance:

An improvement in balance and a reduction in the number of falls cannot be achieved by specific walking training alone. The principle of treatment is a balance and mobility training in the subacute and in the chronic stage after a stroke (ReMoS Working Group, 2015 at 484).

Author's exclusions:

This article is a personal written summary of an official open access publication. This article is not an official document of the guideline working group. It can differ in parts from the original and may contain own interpretations of the author.

Table 7: Improvement of balance (static, dynamic, falls) (ReMoS Working Group, 2015 at 484).

	Subacute	Chronic
B (should):	<ul style="list-style-type: none"> - Intensive gait training without treadmill or - Intensive gait training including treadmill or - Intensive supervised home exercise programme (strengthening, endurance, balance) with progression - Motor relearning programmes 	
0 (can):	<ul style="list-style-type: none"> - Gait training with end-effector or exoskeleton gait trainer, treadmill, if available - Additional ergometer training - Strength/endurance training - Trunk activities on unstable support surface - Additional biofeedback platform with advanced tasks - Acoustic feedback while walking - Needle acupuncture including electroacupuncture during intensive rehabilitation - Additional electroacupuncture - Early orthopaedic shoe 	<ul style="list-style-type: none"> - Gait training with gait trainer or Lokomat if available - Task-related strength training - Exercises on unstable support surface - Exercise programme with systematic reduction of support surface and progression of perturbation - Customised exercise programme (balance, coordination) - Ai Chi (aquatic tai chi) - Movement observation - Additional biofeedback platform with advanced movement components - Additional VR-based training - Combined dual-task-based exercise programme
- B (should not):		<ul style="list-style-type: none"> - Functional electrical stimulation with percutaneous wire electrodes

LITERATURE

[1] **ReMoS Working Group** (Dohle, Quintern, Saal, Stephan, Tholen, Wittenberg). S2e-Leitlinie Rehabilitation der Mobilität nach Schlaganfall (ReMoS). [S2e-Guideline Rehabilitation of Mobility post Stroke (ReMoS)]. *Neurol & Rehabil* 2015, 7: 355-494.

[2] **Moher D, Liberati A, Tetzlaff J, Altman DG**. The PRISMA Group. Preferred Reporting Items for Systematic Reviews and Meta-Analyses: The PRISMA Statement. *PLoS Med* 2009, 6(7).

[3] **Jeremy Howick, Iain Chalmers, Paul Glasziou, Trish Greenhalgh, Carl Heneghan, Alessandro Liberati, Ivan Moschetti, Bob Phillips, and Hazel Thornton**. The 2011 Oxford CEBM Levels of Evidence (Introductory Document).

[4] **Guyatt Gordon H, Oxman Andrew D, Vist Gunn E, Kunz Regina, Falck-Ytter Yngve, Alonso-Coello Pablo et al**. GRADE: an emerging consensus on rating quality of evidence and strength of recommendations *BMJ* 2008, 336:924.

The original full text version and short version of the guideline, written in German, can be downloaded at <http://www.dgnr.de/>.

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