

The compact leg and upper body exerciser with motor

mobi 540



THERA-soft®
connectable + 

THERA-Trainer mobi 540

Item no. A004-419

● = Equipment variant ○ = Accessory

Name	Item no.	Equipment
Base unit mobi	A003-975	●
Foot rests	A003-974	●
Therapy grips	A001-420	●
Control and display unit mobi 2.7" wireless	A006-172	●
Wristbands for therapy grip/safety handle	A002-440	○
Carry bag	A003-976	○
THERA-soft	A006-630	○
Media station incl. PC	A008-119	○
Wireless receiver for PC	A004-020	○

The compact leg and upper body exerciser **with motor**

mobi 540

Product specifications	
Product description	THERA-Trainer mobi is a powered training device using cyclical movement actively and/or passively mobilizing the upper or lower extremities to redevelop muscles, restore motion to joints and improve the cardiovascular function.
Special product features	<ul style="list-style-type: none"> - START/STOP traffic light control: Start and stop training with one press of a button. - Foot rests with safety edge and velcro straps to secure the foot: Ensure a safe workout through optimal support and without foot slippage. - Very good stability and slip resistance: Enables safe training without constant slipping. - Quick and easy cleaning: Due to smooth surfaces, thorough cleaning and disinfection possible. - Tool-free replacement of foot rests and therapy grips: Enables quick and easy adaptation to different treatment settings and therapy goals. - Future-oriented treatment: The THERA-Trainer mobi can be integrated into a THERA-Trainer complete solution at any time.
Place of use	Suitable for use in enclosed spaces in facilities e.g. clinics, nursing homes, therapy practices or at the patient's home
Technical specification	<ul style="list-style-type: none"> - Dimensions (L/W/H): approx. 50 / 46 cm - 56 / 43 cm - Weight: approx. 13.2 kg (with foot rests) or 12.1 kg (with therapy grips) - Electric motor: Nominal power = 100 watts, continuous power = 50 watts - Power connection: 1x 100-240 V, 50/60 Hz - Protection class: II - Degree of protection: Type BF - Protection category: IP 21 - Ambient conditions when in use: Temperature: 5 °C to 40 °C Relative humidity: 15 % to 93 % Rh Air pressure: 700 hPa to 1060 hPa - Sound emission: LpA < 70 dB (A) - Service life: 7 years
Compatibility	THERA-soft
Patient requirements	<ul style="list-style-type: none"> - Patient weight: max. 140 kg - Patient height: 120 - 200 cm
Delivery and commissioning	<ul style="list-style-type: none"> - Packaging (L/W/H): 56 / 46 / 40 cm - Delivery: usually Package via parcel service - Transport and storage conditions: Temperature: -25 °C to 70 °C Relative humidity: 15 % to 93 % Rh Air pressure: 700 hPa to 1060 hPa - Commissioning: by authorised specialists (recommended)
Service / Maintenance	<ul style="list-style-type: none"> - Service: only by authorised service technicians - Maintenance: maintenance-free - Technical inspection: safety inspection every 24 months (manufacturer recommendation)
Registration	<ul style="list-style-type: none"> - Medical device class: IIa - CE marking: CE0297