

Compact active and passive cycling trainer
for upper and lower extremities for home use

mobi



THERA-soft®
connectable + 

THERA-Trainer mobi	
Item no.	U001-600

● = Included ○ = Accessory

Name	Item no.	Included
Base unit	A003-975	●
Foot rests	A003-974	●
Therapy grips	A001-420	●
Control and display unit mobi 2.7" wireless	A006-172	●
Wristbands for therapy grip/safety handle	A002-440	○
Carry bag	A003-976	○
THERA-soft (Annual fee, first year included)	A006-630	○
Media station	A008-119	○
Isolating transformer 115V	A004-002	○

Compact active and passive cycling trainer
for upper and lower extremities for home use

mobi

Product specifications	
Product description	THERA-Trainer mobi is a powered training device using cyclical movement actively and/or passively mobilizing the upper or lower extremities to redevelop muscles, restore motion to joints and improve the cardiovascular function.
Special product features	<ul style="list-style-type: none"> - Foot rests with safety edge and velcro straps to secure the foot: Ensure a safe workout through optimal support and without foot slippage. - Very good stability and slip resistance: Enables safe training without constant slipping. - Quick and easy cleaning: Thorough wipe disinfection is possible thanks to smooth surfaces. - Tool-free replacement of foot rests and therapy grips: Enables quick and easy adaptation to different treatment settings and therapy goals. - Forward-looking treatment: The THERA-Trainer mobi can be integrated into a THERA-Trainer complete solution at any time.
Place of use	Suitable for use in enclosed spaces in facilities e.g. clinics, nursing homes, therapy practices or at the patient's home
Technical specification	<ul style="list-style-type: none"> - Dimensions (LxWxH): approx. 19.7 x 22.1 x 16.9 inch - Weight: approx. 29 lb (with foot rests) or 27 lb (with therapy grips) - Electric motor: Nominal power = 100 watts, continuous power = 50 watts - Power connection: 1x 100-240 V, 50/60 Hz - Protection class: II - Degree of protection: Type BF - Protection category: IP 21 - Ambient conditions when in use: Temperature: 41 °F to 104 °F Relative humidity: 15 % to 93 % Rh Air pressure: 700 hPa to 1060 hPa - Sound emission: LpA < 70 dB (A) - Service life: 7 years
Compatibility	THERA-soft
Patient requirements	<ul style="list-style-type: none"> - Patient weight: max. 308 lb - Patient height: 47.2 - 78.7 inch
Delivery and commissioning	<ul style="list-style-type: none"> - Transport and storage conditions: Temperature: -13 °F to 158 °F Relative humidity: 15 % to 93 % Rh Air pressure: 700 hPa to 1060 hPa - Commissioning: by authorised specialists (recommended)
Service / Maintenance	<ul style="list-style-type: none"> - Service: only by authorised service technicians - Maintenance: maintenance-free - Technical inspection: safety inspection every 24 months (manufacturer recommendation)
Registration	<p>Medical device class: I (USA) CE marking: CE0297 (monitored by certification body)</p>