

Compact active and passive cycling trainer
for upper and lower extremities for home use

mobi



THERA-soft®
connectable + 

● = Included ○ = Accessory

THERA-Trainer mobi

| | |
|----------|----------|
| Item no. | U001-600 |
|----------|----------|

| Name | Item no. | Included |
|--|----------|----------|
| Base unit | A003-975 | ● |
| Foot rests | A003-974 | ● |
| Therapy grips | A001-420 | ● |
| Control and display unit mobi 2.7" wireless | A006-172 | ● |
| Wristbands for therapy grip/safety handle | A002-440 | ○ |
| Carry bag | A003-976 | ○ |
| THERA-soft (Annual fee, first year included) | A006-630 | ○ |
| Media station | A008-119 | ○ |
| Isolating transformer 115V | A004-002 | ○ |

Compact active and passive cycling trainer
for upper and lower extremities for home use

mobi

| Product specifications | |
|----------------------------|--|
| Product description | THERA-Trainer mobi is a powered training device using cyclical movement actively and/or passively mobilizing the upper or lower extremities to redevelop muscles, restore motion to joints and improve the cardiovascular function. |
| Place of application | Suitable for use in enclosed spaces in facilities e.g. clinics, nursing homes, therapy practices or at the patient's home |
| Technical specification | <ul style="list-style-type: none"> - Dimensions (LxWxH): approx. 19.7 x 22.1 x 16.9 inch - Weight: approx. 29 lb (with foot rests) or 27 lb (with therapy grips) - Electric motor: Nominal power = 100 watts, continuous power = 50 watts - Power connection: 1x 100-240 V, 50/60 Hz - Protection class: II - Degree of protection: Type BF - Protection category: IP 21 - Ambient conditions when in use: <ul style="list-style-type: none"> Temperature: 41 °F to 104 °F Relative humidity: 15 % to 93 % Rh Air pressure: 700 hPa to 1060 hPa - Sound emission: LpA < 70 dB (A) - Operational life: 7 years |
| Compatibility | THERA-soft |
| Patient requirements | <ul style="list-style-type: none"> - Patient weight: max. 308 lb - Patient height: 47.2 - 78.7 inch |
| Delivery and commissioning | <ul style="list-style-type: none"> - Transport and storage conditions: <ul style="list-style-type: none"> Temperature: -13 °F to 158 °F Relative humidity: 15 % to 93 % Rh Air pressure: 700 hPa to 1060 hPa - Commissioning: by authorised specialists (recommended) |
| Service and maintenance | <ul style="list-style-type: none"> - Service: only by authorised service technicians - Maintenance: maintenance-free - Technical inspection: safety inspection every 24 months (manufacturer recommendation) |
| Registration | <ul style="list-style-type: none"> - Medical device class: I (USA) - CE marking: CE0297 (monitored by certification body) |