


Cognitive-motor training system

senso



● = Included ○ = Accessory

THERA-Trainer senso 642

Item no.	A008-961
----------	----------

Name	Item no.	Equipment
THERA-Trainer senso incl. - Training platform with handrails - Console with PC and RFID Reader - Tablet - Senso "Play" Software - License period 1 year - HDMI connection cable - Wireless keyboard with mouse - 5 RFID Wristbands	A008-961	●
Media Station	A008-962	○
Senso "Play" Software - License period 1 year (Extension)	A008-964	○
RFID wristband senso	A008-982	○

senso

Product specifications	
Product description	The THERA-Trainer senso is an interactive training and test system for improving cognitive-motor functions. The cognitive-motor training sets new stimuli by interactively coupling with movements. Senso has been proven to improve the performance of elderly and sick people (e. g. reaction and walking speed, balance, stability, concentration, attention, coordination).
Place of application	Suitable for use in enclosed spaces in facilities e.g. clinics, nursing homes, therapy practices or at the patient's home
Technical specification	<ul style="list-style-type: none">- Dimensions (L/W/H): 44.5 / 44.5 / 40.3 inch- Weight: approx. 176.4 lbs- Software information: browser-based, WLAN connection needed- Power connection: 230 V~, 50/60 Hz- Protection class: I- Protection category: IP21- Ambient conditions when in use:<ul style="list-style-type: none">Temperature: 32 °F to 104 °CRelative humidity: 0 % to 60 % Rh- Noise emission: LpA < 60 dB (A)- Operational life: 10 years
Compatibility	Senso "Play" Software
Patient requirements	<ul style="list-style-type: none">- Patient's weight: max. 396.8 lbs
Delivery and commissioning	<ul style="list-style-type: none">- Packaging (L/W/H): 47.2 / 31.5 / 55.1 inch- Delivery: Depending on the equipment, several packages- Transport and storage conditions:<ul style="list-style-type: none">Temperature: 14 °F to 158 °CRelative humidity: 0 % to 60 % Rh- Commissioning: by authorised specialists
Service and maintenance	<ul style="list-style-type: none">- Service: only by authorised service technicians- Maintenance: maintenance-free- Technical inspection: safety inspection every 12 months (manufacturer recommendation)
Registration	<ul style="list-style-type: none">- Medical device class: I- CE marking: CE